

FORMFIT® BACK SUPPORT AIR

The Form Fit Back Support with Air provides compression and support for those suffering from lower back pain. The integrated pneumatic system is intuitive and conveniently located so that adjustments can be made as demands on the body change throughout the day. The double wrap design improves compression, while durable finger loops make applying the brace easy. Made from breathable materials, the Form Fit Back Support provides compression and support for those suffering from lower back pain. The double wrap design improves easy compression, while the durable finger loops make applying the brace easy.

INDICATIONS FOR USE

Lower back conditions that may benefit from gross motion restriction, such as;

- Muscular low back pain
- Lumbosacral strains
- Chronic low back pain
- Sciatica

FEATURES AND BENEFITS

- Self-contained pneumatic system is low-profile and easy to use
- Mesh panels and ventilation holes increase breathability and enhance comfort
- Overlapping double pull strap design provides additional support to the lumbar and abdominal areas
- Flexible stays provide optimum support in the lumbar and abdominal areas
- Gripper strips prevent migration
- Adjustable air bladder provides a custom fit and offers focused compression in the lumbar area

REQUIRED MEASUREMENTS

- Waist circumference

ORDERING INFORMATION

| Size | Waist circumference | Standard |
|----------|---------------------|----------|
| Small | 63 - 76cm | 209133 |
| Medium | 76 - 86cm | 209135 |
| Large | 86 - 97cm | 209137 |
| X-Large | 97 - 107cm | 209138 |
| XX-Large | 107 - 142cm | 209139 |

ACCESSORIES

| Description | Part# |
|-------------|-------|
| Gel pad | GB28 |

