

## INNOVATOR X® POST-OP ELBOW

Benefiting from our proven, patented technologies, Innovator X, with extending struts, can help improve patient outcomes. With superior comfort, clean aesthetics and the most user-friendly hinge in post-op bracing, Innovator X is designed to increase patient satisfaction and compliance.



UPPER EXTREMITY SOLUTIONS

### INDICATIONS FOR USE

- Indications requiring elbow immobilisation and/or range of motion restriction. This may include:
- Stable fractures of the elbow, distal humerus, proximal radius or ulna
- Tendon and ligament injuries or repairs
- Collateral ligament reconstructions
- Elbow hyperextension
- Triceps tendon repair
- Chronic elbow injuries
- Tennis elbow release

### FEATURES AND BENEFITS

- OneTouch™ easy to adjust ROM hinge
- Slide-to-Size™ struts lead to a more customised fit
- Sensil® silicone padding on the strapping system reduces slipping and malleable aluminium cuffs allow for customisation
- Patented Form Fit® 3-dimensional padding and lycra lining provide superior comfort and compliance
- Comfort pad at the neck and shoulder
- Optional anti-supination/pronation arm bar kit is available

### REQUIRED MEASUREMENTS

- Universal, one size

### ORDERING INFORMATION

#### INNOVATOR X

Side	Part#
Left	W-823101
Right	W-823100

### ACCESSORIES

Description	Part#
Innovator X Bar Kit	W-A82310
Replacement softgoods Left	W-823111
Replacement softgoods Right	W-823110